

# Liquid Handling Filter Tips & Boxes Compatability



## Compatability Tip Chart

Please note: This guide is intended for aid in pipette tip selection only and does not form a guarantee.

If you are unsure if the LabCo pipette tips will suit your pipettor, please request your free of charge sample pack.



| Standard Pipette Tips |              |                | Manufacturer Brand |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|-----------------------|--------------|----------------|--------------------|--------|--------|-------|-----------|-----------|--------|-------------------------------|----------|----------|--------|--------|---------|-------|
| Volume Range          | Product Code | Packaging      | LabCo              | Axygen | Biohit | Brand | Eppendorf | Fine Care | Gilson | Labsystems Thermo Finnpipette | Microlit | Nichiryo | Oxford | Rainin | Socorex | Treff |
| 0.1µL to 10µL         | 550.002.500  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                       | 550.002.525  | Racked         |                    | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |
|                       | 550.002.550  | Racked Sterile |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
| 0.5µL to 20µL         | 550.002.505  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                       | 550.002.530  | Racked         |                    | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |
|                       | 550.002.555  | Racked Sterile |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
| 1µL to 200µL          | 550.002.510  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                       | 550.002.535  | Racked         | 3                  | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |
|                       | 550.002.560  | Racked Sterile |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
| 100µL to 1,000µL      | 550.002.515  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                       | 550.002.540  | Racked         | 3                  | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |
|                       | 550.002.565  | Racked Sterile |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
| 100µL to 5mL          | 550.002.520  | Bulk Bag       | 3                  | 3      | 3      |       | 3         |           |        |                               |          |          | 3      |        |         |       |
| 100 µL to 5mL         | 550.002.521  | Bulk Bag       |                    |        |        |       |           |           | 3      |                               |          | 3        |        |        | 3       |       |
| 100µL to 5.5mL        | 550.002.522  | Bulk Bag       |                    |        |        |       | 3         |           |        | 3                             |          |          |        |        |         | 3     |
| 0.5mL to 10mL         | 550.002.523  | Bulk Bag       |                    |        |        |       | 3         |           |        | 3                             |          |          |        |        |         | 3     |

### Tip Ordering:

|                 |              |                | Manufacturer Brand |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|-----------------|--------------|----------------|--------------------|--------|--------|-------|-----------|-----------|--------|-------------------------------|----------|----------|--------|--------|---------|-------|
| Volume Range    | Product Code | Packaging      | LabCo              | Axygen | Biohit | Brand | Eppendorf | Fine Care | Gilson | Labsystems Thermo Finnpipette | Microlit | Nichiryo | Oxford | Rainin | Socorex | Treff |
| 0.1µL to 10µL   | 550.002.600  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                 | 550.002.650  | Racked Sterile |                    | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |
| 0.5µL to 20µL   | 550.002.605  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                 | 550.002.655  | Racked Sterile | 3                  | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |
| 1µL to 200µL    | 550.002.610  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                 | 550.002.660  | Racked Sterile | 3                  | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |
| 100µL to 1000µL | 550.002.615  | Bulk Bag       |                    |        |        |       |           |           |        |                               |          |          |        |        |         |       |
|                 | 550.002.665  |                | 3                  | 3      | 3      | 3     | 3         | 3         | 3      | 3                             | 3        | 3        | 3      | 3      | 3       | 3     |

**Customer Service 1300 522 784**  
[sales@labquipsales.com.au](mailto:sales@labquipsales.com.au)