

Magnetic Stirrer With Digital Display



Round with 20L Capacity

Features

- Stirs up to 20L
- High magnetic field reduces stirrer bar decoupling
- Electronic speed control with constant speed, even during changes in load.

Ordering Information:

Product Code	Description	
400.100.110	Stirrer 20L	Each



Technical Specifications

Product Code	400.100.010
Plate Dimension	135mm Ø
Plate Material	Stainless Steel
Motor Type	DC Brushless Motor
Motor Rating Input	18W
Motor Rating Output	10W
Power	30W
Voltage	240V/60Hz
Stirring Positions	1
Maximum Stirring Volume (H2O)	20L
Maximum Magnetic Stirrer Bar	80mm
Speed Range	100 to 1,500rpm
Speed Display	Digital
Protection Class	IP42
Dimensions	160 x 280 x 85mm
Weight	2.8kg
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%

20 Litre Stirrer with Ceramic Plate

Features

- Stir up to 20L
- Large square, ceramic glass plate
- Resistant to corrosion and easy to clean
- High magnetic field reduces stirrer bar decoupling
- LCD displays both set and actual speeds

Ordering Information:

Product Code	Description	
400.100.020	Stirrer Ceramic Plate 20L	Each



Technical Specifications

Product Code	400.100.020
Plate Dimension	184 x 184mm
Plate Material	Glass Ceramic
Motor Type	DC Brushless Motor
Motor Rating Input	18W
Motor Rating Output	10W
Power	50W
Voltage	240V/60Hz
Stirring Positions	1
Maximum Stirring Volume	20L
Maximum Magnetic Stirrer Bar	80mm
Speed Range	100 to 1,500rpm
Speed Display	Digital, LCD
Protection Class	IP21
Dimensions	215 x 360 x 112mm
Weight	5.3g
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%



Customer Service 1300 522 784
sales@labquipsales.com.au