

Multi-action Mixers with Digital Display



7.5kg Reciprocal Shaker

Features

- 7.5kg Maxi Shaking Weight (including platform)
- Over speed detection and protection
- Durable and Hard Working
- LCD displays speed and timer
- Wide range of platforms for multiple applications



400.102.150
pictured with 400.102.157
Universal Platform & Flasks
all sold separately

Technical Specifications:

Product Code	400.102.150
Voltage	240V/60Hz
Power	30W
Mixing Action	Horizontal
Mixing Diameter	10mm
Max Shaking Weight (with platform)	7.5kg
Motor Type	DC Brushless Motor
Motor Rating Input	28W
Motor Rating Output	15W
Speed Range	100 to 350rpm
Speed Display	Digital, LCD
Timer	Yes
Timer Display	Digital, LCD
Timer Range	1 to 1,199 Minutes
Operation	Continuous / Timed
Dimensions	420 x 370 x 100mm
Weight	13.5kg
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21
RS232	Yes

Ordering Information:

Product Code	Description	Pack Size
400.102.150	Shaker Horizontal 7.5kg	Each

10kg Rocking Mixer

Features

- Smooth see-saw action
- Fixed angle of 9°
- Load up to 10kg
- Digital displays speed and timer
- Wide range of platforms for multiple applications



Technical Specifications:

Product Code	400.105.000
Voltage	240V/60Hz
Power	40W
Motor Type	DC Brushless Motor
Mixing Action	See Saw (Rocking)
Tilt Angle	9°
Maximum Load (including platform)	10kg
Speed Range	10 to 70rpm
Speed Display	Digital, LCD
Timer	Yes
Timer Display	Digital, LCD
Timer Range	1 to 1,199 Minutes
Operation	Continuous / Timed
Dimensions	410 x 360 x 200mm
Weight	9kg
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21

Ordering Information:

Cat. No.	Description	
400.105.000	Shaker See Saw Supplied with 400.105.013 Platform	Each

Customer Service 1300 522 784
sales@labquipsales.com.au