

Orbital Action Mixer 7.5kg with Digital Display



Orbital Action Digital Display Mixer 7.5kg

Features

- 7.5kg Maximum Shaking Weight
- Over speed detection and protection
- DC brushless motor
- LCD displays speed and timer
- Wide range of platforms for multiple applications

Technical Specifications:

Product Code	400.101.150
Voltage	240V/60Hz
Power	30W
Mixing Action	Orbital
Mixing Diameter	10mm
Maximum Shaking Weight (with platform)	7.5kg
Motor Type	DC Brushless Motor
Motor Rating Input	28W
Motor Rating Output	15W
Speed Range	100 to 500rpm
Speed Display	Digital, LCD
Timer	Yes
Timer Display	Digital, LCD
Timer Range	1 to 1,199 Minutes
Operation	Continuous / Timed
Dimensions	420 x 370 x 100mm
Weight	13.5kg
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21
RS232	Yes



400.101.150
Pictured with 400.101.157
Universal Platform & Flasks
all sold separately

Ordering Information:

Product Code	Description	
400.101.150	Shaker Orbital 7.5kg	Each

Racks & Clips for Orbital Motion

Ordering Information: 7.5kg Platforms (320 x 320mm):

Cat. No.	Description	
400.102.155	Universal Platform Horizontally Adjustable Bars	Each
400.102.156	Fixing Clip Platform Fixing Clips Sold Separately	Each
400.102.157	Universal Platform Vertically Adjustable Bars	Each
400.102.158	Dish Platform with Slip Resistant Mat	Each
400.102.159	Replacement Bar For Universal Rack	Each




400.102.113 / 158
Dish Platform



400.102.112 / 157
Universal Rack



400.102.114 / 159
Replacement Bar



400.102.111 / 156
Fixing Clip Platform

Ordering Information:

Cat. No.	Fixing Clip	
400.102.200	25mL Flask	Each
400.102.201	50mL Flask	Each
400.102.202	100mL Flask	Each
400.102.203	200/250mL Flask	Each
400.102.204	500mL Flask	Each



Customer Service 1300 522 784
sales@labquipsales.com.au