

Rotary & Suspension Mixers with Digital Display



Rotary Mixer

Features

- Gentle but effective mix for samples up to 50mL
- Digital control for accuracy
- End over end design
- Wide variety of attachments available

Technical Specifications:

Product Code	400.104.040
Voltage	240V/60Hz
Power	40W
Motor Type	DC Brushless
Tilt Angle	0 - 90°
Speed Range	10 to 70rpm
Speed Display	Digital, LCD
Timer	Yes
Timer Display	Digital, LCD
Timer Range	1 to 1,199 Minutes
Operation	Continuous / Timed
Dimensions	510 x 190 x 260mm
Weight	8kg
Permissible Ambient Temp	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21



400.104.040 Rotisserie Mixer supplied with 400.104.050 (16 x 50mL) Attachment.

Unit has ability to be fully tilted and operated from 0 ≤ 90°.



Ordering Information:

Cat. No.	Description	Each
400.104.040	Rotisserie Mixer supplied with 400.104.050 Platform	Each

Ordering Information - Racks:

Cat. No	Tube Size Vol mL	Each
400.104.045	1.5mL x 48 (Horizontal)	Each
400.104.046	15mL x 24 (Horizontal)	Each
400.104.047	50mL x 24 (Horizontal)	Each
400.104.048	1.5mL x 32 (Vertical)	Each
400.104.049	15mL x 16 (Vertical)	Each
400.104.050	50mL x 16 (Vertical)	Each



Suspension Mixer

Features

- Gentle but effective mix for samples 1.5 up to 50mL
- Digital control for accuracy
- End over end design
- Wide variety of discs available

Ordering Information - Mixer:

Product Code	Description	Pack Size
400.104.015	Disc Mixer supplied with 400.104.022 Platform	Each

Ordering Information - Racks:

Cat. No	Application	Each
400.104.020	For 1.5mL x 60	Each
400.104.021	For 15mL x 16	Each
400.104.022	For 50mL x 8	Each
400.104.029	To hold two discs together	PK4



Technical Specifications:

Product Code	400.104.015
Voltage	240V/60Hz
Power	40W
Motor Type	DC Brushless
Tilt Angle	0 - 90°
Speed Range	10 to 70rpm
Speed Display	Digital, LCD
Timer	Yes
Timer Display	Digital, LCD
Timer Range	1 to 1,199 Minutes
Operation	Continuous / Timed
Dimensions	190 x 190 x 260mm
Weight	8kg
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21

Customer Service 1300 522 784
sales@labquipsales.com.au