

# Vortex & Microplate Mixers



## Vortex Mixer

### Features

- Small, touch or continuous tube mixer
- Vacuum suction feet hold the unit to the bench
- Aluminium cast base avoids unnecessary “walking”
- Versatile unit with a range of tube accessories

### Technical Specifications:

Product Code	400.101.001
Voltage	230V/60Hz
Power	60W
Mixing Action	Orbital
Mixing Diameter	4mm
Motor Rating Input	58W
Motor Rating Output	10W
Speed Range	0 to 2,500rpm
Speed Display	Scale
Operation	Touch / Continuous
Dimensions	127 x 130 x 160mm
Weight	3.5kg
Permissible Ambient Temp	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21



### Ordering Information:

Product Code	Description	Pack Size
400.101.001	Shaker Vortex supplied with 400.101.010 Standard Adapter	Each


## Microplate Mixer

### Features

- Compact and universal
- The Microplate shaker is suitable for 96 well plates or a wide range of tube
- Continuous operation
- DC brushless motor

### Technical Specifications:

Product Code	400.101.050
Voltage	230V/60Hz
Power	20W
Mixing Action	Orbital
Mixing Diameter	4.5mm
Maximum Shaking Weight (with adapter)	0.5kg
Motor Type	DC Brushless Motor
Motor Rating Input	18W
Motor Rating Output	10W
Speed Range	Single Microplate: 0 to 1,500rpm Double Microplate: 0 to 1,000rpm
Speed Display	Scale
Operation	Continuous
Dimensions	260 x 150 x 80mm
Weight	3kg
Permissible Ambient Temperature	5 to 40 ° C
Permissible Relative Humidity	80%
Protection Class	IP21



### Ordering Information:

Product Code	Description	Pack Size
400.101.050	Shaker Microplate supplied with	Each
400.101.061	Microplate (Single) Adapter	

**Customer Service 1300 522 784**  
[sales@labquipsales.com.au](mailto:sales@labquipsales.com.au)