



# The right equipment at the right price!

**Labquip Technologies**

## Fan Forced Drying Oven with LCD Display

Forced-air circulation combined with digital temperature control achieves dependable and uniform environment for a variety of laboratory applications.

Corrosion resistant stainless steel interior finish with toughened double glazing observation window.

Gentle downward flow along chamber wall provides less disruption to sample material.

Full 12 month warranty.



**Built-in temperature calibration**

**User friendly soft-touch control panel**

**Auto overheat protection**

**CE compliant**

Model	F0-19040	F0-19070	F0-19140	F0-19240	F0-19440	F0-19640
Power	220V/50Hz single-phase			220V/50Hz single-phase		480V/50Hz three-phase
Output	1000 watts	1500 watts	2000 watts	2500 watts	3500 watts	6400 watts
Alarm Range	±5°C					
Uniformity	±1°C					
Temperature Range	Ambiant - 250°C					
Chamber Size (mm)	356x300x356	460x356x460	560x406x560	610x510x762	686x660x889	787x787x991
Machine Size (mm)	508x535x661	610x584x762	686x686x864	762x686x1067	838x813x1270	965x991x1321
Chamber Capacity (m <sup>3</sup> )	0.04	0.07	0.14	0.23	0.41	0.62

All equipment sold is covered by statutory 12 month warranty for performance

[www.consumer.vic.gov.au](http://www.consumer.vic.gov.au)

DRIVEN BY  
**Labquip Technologies**

[sales@labquipsales.com.au](mailto:sales@labquipsales.com.au)  
**Tel: 1300 LABQUIP (1300 522 784)**  
visit: [labquip.com.au](http://labquip.com.au)