



The right equipment at the right price!

Labquip Technologies

Digital Forced Air Incubator Oven

Uses microcontroller to enhance temperature precision, accuracy and uniformity.

The memory function saves you time by retaining all of your settings when the power is down and will resume previous oven operation once powered on again.

1 - 9999 min timer function.



Built-in temperature calibration

Temperature out-of-range alarm

Auto overheat protection

CE compliant

Model	FO-19013	FO-19023	FO-19053
Output	500 watts	500 watts	800 watts
Power	220V/50Hz single phase		
Temperature range	Ambient °C to 200°C		
Plate temp. uniformity	±1°C		
Alarm range	±5°C		
Chamber Size (mm)	250x230x250	300x280x280	450x290x460
Machine Size (mm)	533x432x432	610x508x457	711x559x533
Chamber capacity (m ³)	0.014	0.025	0.054

All equipment sold is covered by statutory 12 month warranty for performance

www.consumer.vic.gov.au

Labquip DRIVEN BY
Technologies

sales@labquipsales.com.au
Tel: 1300 LABQUIP (1300 522 784)
visit: labquip.com.au