



The right equipment at the right price!

Labquip Technologies

Thermostatic Water Baths

Range of laboratory thermostatic water baths.

From single hole to 8 holes across two rows.

Stainless steel inner tank.

Touch button setting control, digital display with fixed temperature settings, timer and auto stop.

Deviation correction, menu key lock, power failure backup, setting memory on shut down, PID adjustment and overheat alarm.



| Model | WB-3-1 1 hole | WB-7-2 1 row, 2 holes | WB-13-4 1 row, 4 holes | WB-19-6 1 row, 6 holes | WB-12-4 2 rows, , 4 holes | WB-18-6 2 rows, 6 holes | WB-20-8 2 rows, 8 holes |
|--------------------------|---------------------------------------|--------------------------|---------------------------|---------------------------|------------------------------|----------------------------|----------------------------|
| Heating mode | Bottom heating natural conection | | | | | | |
| Temp. range | R1+5 ~ 100°C | | | | | | |
| Temp resoluton range | 0.1°C | | | | | | |
| Temp motion | ±0.5°C | | | | | | |
| Temp uniformity | ±1°C | | | | | | |
| Power | 300W | 600W | 1000W | 1500W | 800W | 1200W | 2000W |
| Timer | 0 ~ 999 minute (with timing function) | | | | | | |
| Inner chamber size (mm) | 160x160x130 | 305x160x130 | 600x165x130 | 900x305x130 | 305x305x130 | 470x305x130 | 600x290x130 |
| Exterior Size (mm) | 200x210x210 | 340x210x210 | 635x165x210 | 930x210x210 | 340x350x210 | 510x350x210 | 635x340x210 |
| Packing Size (mm) | 270x270x275 | 405x275x275 | 715x260x275 | 1000x260x275 | 415x415x275 | 575x410x275 | 700x415x275 |
| Internal volume | 3.3L | 6.3L | 12.8L | 19.3L | 12L | 18.6L | 20.6L |
| Cover aperture | 120mm | | | | | | |
| Nett weight/Gross weight | 3/3.5kg | 4/5kg | 7/8kg | 8/10kg | 6/7kg | 8/9kg | 9/11kg |

All equipment sold is covered by statutory 12 month warranty for performance

www.consumer.vic.gov.au

Labquip Technologies

sales@labquipsales.com.au
Tel: 1300 LABQUIP (1300 522 784)
visit: labquip.com.au