



The right equipment at the right price!

LABQUIP Sciences

Ultra Sonic Cleaner

30L Digital Ultrasonic Cleaner

This commercial-grade ultrasonic cleaner is built with sturdy 304 stainless steel, making it an efficient machine for daily cleaning work. Our ultrasound jewellery cleaner machine will emit 40kHz ultrasonic waves to deeply wash dirt on your valuables' surfaces and crevices.

With heating and timing functions, important items can be cleaned thoroughly in minutes. LabSonic ultrasonic cleaners are ideal and suitable for cleaning laboratory glassware, fine metal objects and jewellery.

Full 12 month warranty.




Built in temperature calibration

User friendly soft-touch control Panel

Auto overheat protection

CE compliant



Adjustable Time/Temperature

Time is adjustable from 0 to 30 minutes and temperature can be set within RT-176°/RT-80° displayed on LED screens. The heating function dramatically improves cleaning efficiency

Model	LS-30
Power	220V/50Hz
Output	40kHz
Uniformity	±1°C
Temperature Range	Ambiant - 80°C
Bath Size (mm) ID	475 x 272 x 195
Unit Size (mm) OD	517 x 312 x 302
Volume (ltr)	30

Cat No	Size	Vol	\$
LS-30	470x270x150	30ltr	2700.00
LS-15	350x320x150	15ltr	1900.00
LS-10	320x260x150	10ltr	1150.00
LS-6	320x170x150	6ltr	750.00

LABQUIP Sciences

sales@labquipsales.com.au
Telephone 1300 522 784