

Chemicals and Life Science Scientific Laboratory Equipment

LABQUIP Sciences

20lt High Torque Over-Head Mixer

Included: Motor stand
mixing stirrer blade



Robust and reliable suitable for Emulsions, Creams,
Mixing Liquids and Paints/Dyes

LCD display for precise monitoring of set and actual
speeds, a wide range of speeds.

Count down Timer.

Speed Range 50 to 2200rpm with control accuracy
of ± 3 rpm.

Brushless DC motor for long life, maintenance free.

Smooth operation prevents accidental spillage
and splashing.

Self Adjusting speed even with changes in
viscosities of the samples.

Auto-restart function in case of power loss.

Safety circuits allow for safe stop function in
anti-stall or overload conditions.

Specifications	OLS20-Pro
Max stirring quantity (H ₂ O)	20L
Power	70W
Motor type	Brushless DC motor
Motor input power	60W
Motor output power	50W
Speed range	50-2200rpm
Speed display accuracy	± 3 rpm
Speed display	LCD
Torque display	LCD
Voltage	100-220V, 50/60Hz
Overhead protection display	LED lights
Max. torque	40Ncm
Viscosity max.	10000.Pas
Chuck range diameter	0.5 - 10mm
Dimensions (WxDxH)	83x220x186mm
Weight	2.8kg
Protection class	IP21
Permissible ambient temp and humidity	5-40°C, 80%RH
Data connector	RS232

LABQUIP Sciences

sales@labquipsales.com.au
Tel: 1300 LABQUIP (1300 522 784)
visit: labquip.com.au