

The right equipment at the right price!

LABQUIP Sciences

Benchtop pH Meter

Benchtop laboratory pH Meter mV Tester. Temperature compensation with digital screen LCD. Model PHS-25

Supplied with:
pH Meter
Electrode
Power Supply
Stand



LCD digital screen

Benchtop pH meter



Specification	Value
Measuring Range	pH0.00-14.00pH; mV0-1800mV
Resolution	0.01pH, 1mV, 0.1°C
Manual Temp. Compensation	0-60°C
Basic Error of Electronic Unit	pH0.01pH; mV0.1% (FS)
Instrument Basic Error	0.02pH
Electronic Repeatability Error	pH0.01pH; mV1mV
Instrument Repeatability Error	≤0.01pH
Electronic Stability	0.01pH/3h
Dimensions	220x200x85mm (lxbxh)

LABQUIP Sciences

sales@labquipsales.com.au
Telephone 1300 522 784