



# The right equipment at the right price!

**LABQUIP** Sciences

## Ohaus Precision Balance



### Ohaus Precision Balance (SKX)

11 models available, capacities from 120g x 0.001g to 8200g x 1g - round or square stainless steel weighing platforms.

The Scout General SKX is complete with a large LCD display with backlight and 4 button control enabling efficient operation and ease of use.

Overload protection rated at 10 times the capacity, making it the perfect partner for analytical, research and teaching applications.

Model	Spec	\$
SKX123	120g x 0.001g	780.00
SKX222	220g x 0.01g	590.00
SKX422	420g x 0.01g	790.00
SKX622	620g x 0.01g	860.00
SKX1202	1200g x 0.01g	950.00
SKX2202	2200g x 0.01g	1020.00
SKX421	420g x 0.1g	500.00
SKX621	620g x 0.1g	590.00
SKX2201	2200g x 0.1g	660.00
SKX6201	6200g x 0.1g	770.00
SKX8200	820g x 1g	590.00

**Large LCD digital display**

**Simple setup**

**Portable**

Removable stainless steel weigh pan

Sealed front panel

Stability indicator

Overload/underload Protection

Low battery indicator

Auto shut off

AC or 4 x AA battery operation

RS232, USB host, USB device or ethernet - optional

1 year manufacturers warranty

No	Option	\$
30045641	SF40A Impact printer	760.00
30268987	Stacking & storage cover set of 6*	109.00
30268988	Stacking & storage cover set of 1*	30.00
30269021	Carrying case	205.00
30268986	Ethernet interface kit	325.00
30268984	USB interface kit	174.00

\*Not suitable for SKX123

**LABQUIP** Sciences

sales@labquipsales.com.au  
Telephone 1300 522 784